

General TENS Information

Electrode Placement in TENS therapy:

- Place directly over the area from which the pain is originating
- “Encircle” the area of pain
- On specific points such as trigger or acupuncture points
- Proper instruction on electrode placement will be given by a Golden State Medical representative at the time the TENS unit is issued.

Should you experience a problem with your TENS device:

- Check whether the Amplitude control is set to the appropriate form of therapy and adjust accordingly
- Check whether the wires are correctly connected to the device and are inserted completely into the sockets
- Check whether the power indicator light is illuminated. If necessary, insert a new battery
- DO NOT attempt to repair a defective device! Contact Golden State Medical for replacement/repair procedures