

## **GENERAL TENS INFORMATION**

### **Electrode Placement in TENS therapy:**

- Place directly over the area from which the pain is originating.
- “Encircle” the area of pain.
- On specific points, such as trigger or acupuncture points.
- Proper instruction on electrode placement will be given by a Golden State Medical, Inc. representative at the time the TENS unit is issued.

### **Should you experience a problem with your TENS unit:**

- Check to see that the Amplitude Control is set to the appropriate form of therapy and adjust accordingly.
- Check to see if the wires are correctly connected to the unit and are inserted completely into the sockets.
- Check to see that the power indicator light is illuminated. If necessary, insert new batteries.
- DO NOT attempt to repair a defective unit! Contact Golden State Medical, Inc. for replacement/repair procedures.