

TREATMENT KITS



ionto+*HI-PER*
plus

Iontophoresis Electrodes *plus* Application kit

HI-PERFORMANCE
for Higher mA protocols



Active pad
Drug reservoir electrode
(placed over treatment area)



Dispersive pad
Return electrode
(completes circuit)



- Greater adhesive area eliminates pad movement
- Unique Silver-Silver Chloride conductor (superior to carbon)
Higher Conduction=More Drug Delivery
- Buffering agents for additional pH stabilization

- Synthetic Gel offers better adhesion than Karaya.
Provides maximum current dispersion
Increases patient comfort level
- Synthetic Buffering provides consistent pH balancing
Helps prevent the effects of hydrolysis
- No water needed to activate the adhesive
No chance of contamination by water

TREATMENT KITS

SMALL
(1.5cc/3.0mA
2.75"x2.75"

BUTTERFLY



MEDIUM
(2.5cc/4.0mA
3.25"x3.25"



LARGE
(4.0cc/4.0mA
3.5"x3.5"



(2.0cc/3.0mA
3.375"x3.75"

Part #	Description
201-452	Butterfly (12 Treatments)
201-455	Small (12 Treatments)
201-449	Medium (12 Treatments)
201-453	Large (12 Treatments)